



✓ = Available
 For 1-piece abutments, either brand screws or CMC screws may be used. Implant brand screws are not provided by CMC. All registered marks are owned by the respective manufacturers and are not the property of Custom Milling Center.



Implant Brand	Platform Type	Titanium	Hybrid	Screw Torque (N-CM)
Astra Tech	OsseoSpeed 3.5/4.0	✓	✓	25
	OsseoSpeed 4.5/5.0	✓	✓	25
Biomet 3i™	Certain 3.4	✓	✓	20
	Certain 4.1	✓	✓	20
	Certain 5.0	✓	✓	20
	Certain 6.0	✓	✓	20
Camlog*	Φ3.3	✓	✓	20
	Φ3.8	✓	✓	20
	Φ4.3	✓	✓	20
	Φ5.0	✓	✓	20
	Φ6.0	✓	✓	20
Conelog*	Φ3.3		✓	20
	Φ3.8		✓	20
	Φ4.3		✓	20
	Φ5.0		✓	20
Friadent XiVE™	Frialit 3.0			24
	Frialit 3.4		✓	24
	Frialit 3.8		✓	24
	Frialit 4.5		✓	24
	Frialit 5.5		✓	24
	Frialit 6.5		✓	24
Intra-Lock*	SQ (Straight 3.3/Taper 3.75)			30
	Standard (Straight 4.0/Conic 4.3)		✓	30
	Wide (Straight 4.75/Conic 4.3, 5.0)		✓	30
Implant One Platform	Series 100	✓	✓	20
	Series 200	✓	✓	20

LifeCore	Internal Hex SDI	✓	✓	30
	Internal Hex RDI 4.5/4.75	✓	✓	30
Megagen *Non Engaging Options Available	AnyRidge Internal *		✓	30
	AnyOne Internal *		✓	30
	EZ Plus Internal Small		✓	30
	EZ Plus Internal Reg.		✓	30
	ExFeel Internal		✓	30
	Rescue Internal		✓	30
Neodent	CM Exact	✓	✓	15
Neoss*	3.25	✓	✓	32
	3.5-5.5	✓	✓	32
Nobel Biocare®	Branemark RP	✓		35
	Replace NP	✓	✓	35
	Replace RP	✓	✓	35
	Replace WP	✓	✓	35
	Replace 6.0	✓	✓	35
	Active NP	✓	✓	35
	Active RP	✓	✓	35
Osstem*	GS Mini (Φ3.5)		✓	30
	GS Standard (Φ4.0/Φ4.5/Φ5.0)		✓	30
	SS Regular (Φ4.8)		✓	30
	SS Wide (Φ6.0)		✓	30
Straumann*	synOcta™ ITI NN	✓		35
	synOcta™ ITI RN	✓	✓	35
	synOcta™ ITI WN	✓	✓	35
	Bonelevel NC (Φ3.3)	✓	✓	35
	Bonelevel RC (Φ4.1/Φ4.8)	✓	✓	35
Sweden & Martina	Premium 3.3	✓	✓	20
	Premium 3.8	✓	✓	20
	Premium 4.25	✓	✓	20
	Premium 5.0	✓	✓	20
	Shelta 3.8	✓	✓	20
	Shelta 4.25	✓	✓	20
	Shelta 5.0	✓	✓	20
Zimmer	Tapered Screw-vent 3.5	✓	✓	30
	Tapered Screw-vent 4.5	✓	✓	30
	Tapered Screw-vent 5.7	✓	✓	30